



TSCD[®]-Q

STERILE TUBING WELDER

TERUMOBCT
Unlocking the Potential of Blood

DISCOVER POLYURETHANE TUBING
CONNECTIONS BASED ON PROVEN
TERUMO BCT TECHNOLOGY



SIMPLE

STERILE CONNECTIONS IN A FEW EASY STEPS



CONSISTENT

STRONG, RELIABLE WELDS, EACH AND EVERY TIME



POLYURETHANE-SPECIFIC

VALIDATED FOR POLYURETHANE
TUBING CONNECTIONS

Learn more at [TERUMOBCT.COM/TSCDQ](https://www.terumobct.com/tscdq)



THE TSCD-Q STERILE TUBING WELDER FROM TERUMO BCT

POLYURETHANE TUBING CONNECTIONS YOU CAN TRUST

Built on a legacy of proven thermal welders and cell processing devices from Terumo BCT, the TSCD-Q sterile tubing welder creates strong, sterile polyurethane tubing connections.



A STERILE PROCESS TO MAINTAIN A FUNCTIONALLY CLOSED SYSTEM

This compact, portable welder allows the attachment and removal of components through a sterile process, maintaining a functionally closed system and providing protection from outside contamination.



THE FLEXIBILITY TO WELD WET-TO-DRY, WET-TO-WET AND DRY-TO-DRY CONNECTIONS

Designed for use with 0.187 in/4.75 mm O.D. (outside diameter) tubing with a wall thickness of 0.031 in/0.787 mm. The TSCD-Q gives you the capability to create wet-to-dry, wet-to-wet and dry-to-dry connections.

**Learn more at TerumoBCT.com/TSCDQ
Schedule a demo today at cellprocessing@terumobct.com**

**Not available in all countries.*

TERUMOBCT

Terumo BCT is a global leader in blood component, therapeutic apheresis and cellular technologies. We believe in the potential of blood. This belief inspires our innovation and strengthens our collaboration with customers.

UNLOCKING THE POTENTIAL OF BLOOD | TERUMOBCT.COM

Terumo BCT, Inc.

10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA

USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT (Asia Pacific) Ltd.

Room 3903-3903A, 39/F
ACE Tower, Windsor House
311 Gloucester Road
Causeway Bay, Hong Kong

Phone: +852.2283.0700
Fax: +852.2576.1311

Terumo BCT Japan, Inc.

20-14, 3-chome
Higashi Gotanda, Shinagawa-ku
Tokyo, 141-0022
Japan

Phone: +81.3.6743.7890
Fax: +81.3.6743.9800

Terumo BCT Europe N.V.

Europe, Middle East and Africa
Ikaroslaan 41
1930 Zaventem
Belgium

Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

Terumo BCT Latin America S.A.

La Pampa 1517–12th Floor
C1428DZE
Buenos Aires
Argentina

Phone: +54.11.5530.5200
Fax: +54.11.5530.5201